

Patient Specific Disposable Slings & Belts



Medcare's Patient Specific Disposable Slings and Belts are designed for use by a single patient and help reduce the risk of cross-infection and -contamination. Although manufactured to be disposable, they can be used multiple times on the same patient. Once washed, patient specific slings are no longer safe to use with patients and the wash indicator (shown below) will be visible. If you see the wash indicator on a sling this signifies it is time to dispose of the sling. **Medcare slings are the first and only slings manufactured with Microban® antimicrobial product protection.**



Medcare's patient specific slings can be used with most mobile ceiling lift styles and fit most 2, 4 and 6 point spreader bars.

Sling & Belt Sizing

Use this sizing chart as a general guide. Keep in mind that two patients that are the same weight may have a completely different body shape and size and may require different size slings or belts.

Sling Sizes

Weight	Sling Size	Strap Color
>100 lbs.	Small	Black
100-210 lbs.	Medium	Blue
210-310 lbs.	Large	Yellow
310-450 lbs.	X-Large	White
450-1000 lbs.	Bariatric	Purple

Belt Sizes

Belt Size	Belt Waist Size	Strap Color
Small	24" to 48"	Black
Medium	30" to 54"	Blue
Large	36" to 60"	Yellow
X-Large	42" to 66"	White

Medcare Slings Contain Antimicrobial Product Protection

Microban antimicrobial protection works 24/7 for an added level of cleanliness protection. Microban inhibits the growth of microbes such as bacteria, mold and mildew that can cause stains, odors and product deterioration on our slings.




Features

- Slings hold up to 1000 lbs.
- 100% breathable and may be left under patient
- May be used for bathing, showering or in the pool as slings will not rip or fall apart when wet
- Can be spot-cleaned or spray disinfected
- Agitation when washing will activate wash indicator signifying time for disposal



Activated wash indicator



MADE IN THE USA



Patient Specific Disposable Slings & Belts (cont.)

Med
HANDLE WITH CARE.

Disposable Care Sling

The Care Sling is our most popular sling style and is designed to fit the needs of most patients.

Item Numbers:

- 600020 – Medium
- 600022 – Large
- 600017 – X-Large
- 600015 – Bariatric Sling



Disposable Care Sling with Head Support

This sling is ideal for patients with weak or poor head control.

Item Numbers:

- 600030 – Medium Sling
- 600032 – Large Sling
- 600034 – X-Large Sling
- 600036 – Bariatric Sling



Disposable Commode Sling

Commode slings are designed to provide easy transfers to and from commodes or toilets and for hygienic purposes. An enlarged opening makes utilizing the commode quicker and easier than standard slings, with minimal interference of clothing and sling material. An adjustable belt attaches around the waist for added security and comfort during transfers.

Item Numbers:

- 600071 – Small Commode Sling
- 600072 – Medium Commode Sling
- 600073 – Large Commode Sling
- 600074 – X-Large Commode Sling



Disposable Repositioning Sling

The Repositioning Sling makes turning or boosting a patient very easy. It is breathable and meant to be left on the bed under the sheet. Dimensions: 55.5" w x 80.5" l

Item Number:

- 600210 – Repositioning Sling



Disposable Repositioning Sling

Make sure you are using slings & belts with the highest level of protection. Look for this tag.

Disposable Stretcher & Morgue Sling

For patients that must be kept in a supine position, the Stretcher Sling allows for easy transfer of the patient from bed to gurney or operating table. This sling is breathable and may be left under a patient. For morgue use, it eliminates manual lifting of cadaver and enables easy transfer between storage units and autopsy table. Dimensions: 47.5" w x 80.5" l

Item Number:

- 600200 – Stretcher Sling
- 600220 – Morgue Sling



Extended Hanger Bar

The extended hanger bar is for use with Morgue and Stretcher slings and is rated at 1000 lbs.

Item Number:

- 400609 – For use with 1" strap Medicare Ceiling Lift
- 400604 – For use with Medicare Care Lift
- 400603 – For use with 1 3/4" strap Medicare Ceiling Lift



Disposable Limb Sling

Limb Slings are designed to assist in lifting a leg or arm for cleaning, changing of wound dressings, etc. They are available in two widths, 6" and 12".

Item Number:

- 600250 – 6" Limb Sling
- 600275 – 12" Limb Sling



Disposable Care Belt

This belt is for weight-bearing patients who have more consistent and reliable motor functions.

Item Numbers:

- 600080 – Medium Belt
- 600082 – Large Belt
- 600084 – X-Large Belt



MEDCARE
WITH ANTIMICROBIAL PROTECTION

THIS SLING PROTECTED BY

Microban
antimicrobial product protection

Cleaner. Fresher. All the Time.™



Medcare Products

151 East Cliff Road
Burnsville, MN 55337
Toll-free: 800-695-4479
Fax: 952-894-7153
www.medcarelifts.com

© 2009 Medicare Products Inc.
Specifications are subject to change without notice. All products referred herein are registered trademarks or trademarks of their respective holders. MICROBAN is a registered trademark of Microban Products Company.