



# NOR-LAKE SCIENTIFIC

## GENERAL PURPOSE LABORATORY FREEZER

Designed to meet the demanding requirements of scientific and laboratory applications. Industrial quality cabinet with a heavy duty refrigeration system to maintain your valuable products.


### Standard Product Features:

- Adjustable control range from -10°C to -25°C  
Factory preset at -20°C.
- Internal chamber thermometer
- The exterior front, sides, back and door are white baked on enamel
- The exterior top and bottom are galvanized steel
- Interior liner and interior door panel are smooth aluminum
- Cabinet and doors are foamed-in-place with CFC Free high density polyurethane foam insulation
- Heavy duty spring loaded hinges
- Key lock door
- Full length stainless steel pull handles
- Four casters (2 locking)
- 3 epoxy coated wire shelves
- Manual defrost
- Complete top mount refrigeration system
- Cold wall evaporator
- Air-cooled condensing unit
- Self closing solid door
- Magnetic door gasket
- Door opening heaters
- Warranties: 18 month parts and labor, 5 year compressor (US and Canada) 18 months parts (International)
- UL and C-UL listed; CEC California Energy Commission; Energy Star



### Optional Features:

- Legs in lieu of casters
- Extra shelves
- Stainless steel exterior (with galv. back)
- Stainless steel interior
- Temperature chart recorder
- Chart paper
- Digital thermometer monitor with alarm (field installed)
- Stainless steel drawers (8 per door)



**NOR-LAKE SCIENTIFIC**  
A Division of  
**NOR-LAKE, INCORPORATED**  
727 Second Street  
P.O. Box 248  
Hudson, Wisconsin 54016  
800-477-5253  
715-386-2323  
**800-388-5253 PARTS**  
715-386-4290 FAX  
[www.norlake.com](http://www.norlake.com)